

MESOTHERAPY/ HYPODERMIC NEEDLE

For hyaluronic acid ,botulax,
red ampoule solution ,hyaluronidase



Catalog

- 01** Product Introduction
- 02** Product Specifications
- 03** Products Advantage
- 04** Products Application
- 05** Products Packaging

01 Product Introduction

INTRODUCTION

1. Mesotherapy uses injections to deliver compounds which may include hyaluronic acid, growth factors and vitamins /minerals directly into the skin to nourish and improve its appearance.
2. This treatment has a widerange of applications including collagen stimulation, hair follicle stimulation, cellulite reduction and improvement of blood circulation.
3. General injection needle is used for subcutaneous, muscular, intravenous or liquid injection in human body





02 Product Specifications


SPECIFICATIONS

	GAUGE	NEEDLE LENGTH	COLOR
	18G	30mm	Pink
	25G	13mm	Orange
	27G	13mm	Gray
	30G	4mm	Yellow
		13mm	
		25mm	
	31G	13mm	Purple
	32G	4mm	Green
		13mm	
	34G	4mm	Orange

02 Product Specifications

SPECIFICATIONS

Mesotherapy/Hypodermic Needle
color example:



03 Products Advantage

ADVANTAGE

- ✓ Medical stainless steel needle Sleek and sharp
- ✓ Sharp needle tip, no blockage,
- ✓ less pain, less tissue damage
- ✓ Compatible with different types of syringes





04 Products Application

APPLICATION

Suitable for all parts of the body, can achieve the best results.

Medical grade

Medical



Fat dissolving
Face lift

Face lift

Water light whitening

Whitening



Mesotherapy Needle

NEEDLE SPEC.

30G	30G	30G
30x10mm	30x15mm	30x20mm
30G	28G	27G
30x10mm	30x15mm	30x20mm
25G	22G	21G
30x10mm	30x15mm	30x20mm
18G	24G	DEM
30x10mm	30x15mm	

Caution:

- Sterile, for single use only.
- Sterility is guaranteed unless the individual package is opened or damaged.
- Do not expose to direct sunlight. Store in a cool & dry place.
- Do not store at extreme temperatures or high humidities.

Caution:

- Do not store at extreme temperatures or high humidities.
- Do not expose to direct sunlight. Store in a cool & dry place.
- Individual package is opened or damaged.
- Sterility is guaranteed unless the
- Sterile, for single use only.

30G

100pcs

DISPOSABLE NEEDLE

Through A Special Process Better
Toughness ,Cleaner Inner Wall
Smoother Indented Groove

CE 0197
STERILE EO

SIZE: 30Gx38mm
LOT NO.: 200200008
MFG DATE: 200200008
EXP DATE: 200700007



05 Products Packing

Korean brand Packaging:



Custom packaging:



Because of the development of the company,
our business further expands

We sincerely recruit

INVESTMENT

wholesalers / The agent

Welcome to join us~

